

KM-5000TP Series

Reclining, lightweight aluminum detachable-arm wheelchair/transport



- Folding 6061 T-6 Aircraft-grade aluminum frame in black powder coat finish
- Medicare Code: K0004+E1226+E0966
- Backrest Reclines 90°~163° and includes anti-tippers
- Full length padded detachable armrests with side panels
- Swing-away detachable elevating footrests
- Adjustable length leg supports & footplates
- Adjustable-height head pillow
- 14” Mag wheels with flat free polyurethane tires
- 7” x 1” Flat free front polyurethane casters
- Adjustable toggle wheel locks
- High-grade flame retardant, breathable nylon upholstery
- Limited lifetime warranty on frames and cross-braces, one year on parts
- One of only two lightweight recliners available in the USA
- Also Available in Large Manual Rear Wheel option, (KM-5000F), please call for more info
- Frame Color: Black

(*Dimensions shown are for Small 14” Wheel Model)

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM-5000TP-16	16	17	20	33 (19+15)	51	25	33	250
KM-5000TP-18	18	18	20	33 (19+15)	51	27	33	250
KM-5000TP-20	20	18	20	33 (19+15)	51	29	34	300
KM-5000TP-22	22	18	20	33 (19+15)	51	31	34	350

Note: Specifications are subject to change without notice



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335