

# KM-8020F-WB

*Aluminum high strength, lightweight  
W/hand brakes & flip-back-arm  
wheelchair*



- Folding 6061 T-6 Aircraft-grade aluminum frame with powder coat finish
- Medicare: K4
- Padded flip-back full sized padded armrests with side panels
- Hand wire brakes for companion and user
- Adjustable rear axle position
- Swing-way detachable footrests with Velcro mount leg strap
- Adjustable length anti-slip footplates
- 24" x 1" high profile flat free polyurethane tires
- High-grade flame retardant, breathable nylon upholstery

Model No.	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
KM 8020F-WB-18	18	16	20	16	38	26	29	250
KM-8020F-WB-16	16	16	20	16	38	24	29	250

Model No.	Description
KM-8020F-WB-18	18", aluminum, high strength/lightweight, hand brakes, full size flip-back-arms, swing away footrests
KM-8020F-WB-18-E	18", aluminum, high strength/lightweight, hand brakes, full size flip-back-arms, elevating footrests
KM-8020F-WB-16	16", aluminum, high strength/lightweight, hand brakes, full size flip-back-arms, swing away footrests
KM-8020F-WB-16-E	16", aluminum, high strength/lightweight, hand brakes, full size flip-back-arms, elevating footrests

Frame Color: **Black**

Common Options: \*Not available on all models \*\* Price & specifications are subject to change without notice.

- Elevating Leg Rests • Anti-Tippers • Brake Extension • IV Pod Attachment • Oxygen Tank Holder • Seat Belt • Heel Loops • Leg Straps • Quick Release Axle • Spoke Wheels • Pneumatic Tires • Run Flat Inserts • Desk Tray • Carrying Pouches •



**KARMAN HEALTHCARE Inc.**

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order.1-800-805-2762 Tel.626-581-2235 Fax.626-581-2335