

LT-1000HB

Aluminum high strength, ultra lightweight hand brake companion chairs



- Durable & lightweight aluminum frame
- Padded fixed full length armrests
- Folding backrest for easy storage
- 19" seat with carry loops for transporting convenience
- Seat belt comes standard w/ transport!
- Swing-away footrests w/ heel loop
- 8" front caster, 12 ½" rear solid polyurethane tire
- Hand loop brake system for attendant
- Manual lock brakes
- High-grade flame retardant, breathable nylon upholstery
- Blue or Burgundy frame available

	seat width (in.)	seat depth (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall Open width (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
LT-1000HB	19	16	19	19	40	21.5	23	250

Model No.	Description
LT-1000HB	19" seat width, folding backrest, fixed full arms, loop brakes, heel loops, and swing-away footrests
LT-1000HB-E	19" seat width, folding backrest, fixed full arms, loop brakes, heel loops, and elevating leg rests

** Price & specifications are subject to change without notice.

OPTIONS:

- Elevating Leg-Rests • Leg Straps



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry. CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335