

LT-700T

Standard lightweight detachable arm wheelchairs



- Medicare code K0003
- Detachable/Padded desk length armrest with alloy side guard
- Swing-away/Length adjustable footrest with heel loop footplate
- Push to lock wheel brakes
- High-grade/Padded flame retardant nylon upholstery
- Dual-Axle/Seat height adjustable
- 24" x 1" Mag wheels with polyurethane tires
- Extra-durable forks
- 8" x 1" Caster wheel
- Carry pocket on backrest
- Seat width: 16", 18"
- Frame color: Silver vein
- Weight capacity: 250 lbs

SPECIFICATIONS (subject to change without notice) :

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
LT-700NT	16	16	10	17 or 19	17	36	23	11	31	36	250
LT-700T	18	16	10	17 or 19	17	36	25	11	31	36	250

OPTIONIONAL ACCESSORIES:

- Elevating Legrests - Anti-Tippers - Brake Extension - IV Pod Attachment - Oxygen Tank Holder - Seat Belt - Leg Straps
- Headrest - Carrying Pouches



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue . City of Industry . CA 91748-1418

Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335