



**KARMA HEALTHCARE**  
 19255 San Jose Ave.  
 City Of Industry, CA 91748  
**Tel:** 626-581-2235 **Fax:** 626-581-2335  
**Toll Free:** 800-805-2762



QR Reader Smartphone Info



**S-ERGO ATX SERIES ORDER FORM**

Submitting for  QUOTE  ORDER

**Mark For:** \_\_\_\_\_  
**Date:** \_\_\_\_\_  
**Dealer Acct#:** \_\_\_\_\_  
**Dealer:** \_\_\_\_\_  
**Dealer Contact:** \_\_\_\_\_  
**Dealer Address:** \_\_\_\_\_  
**Dealer City:** \_\_\_\_\_ **ST:** \_\_\_\_\_ **ZIP:** \_\_\_\_\_  
**Dealer Phone:** (\_\_\_\_) \_\_\_\_\_ **FAX:**(\_\_\_\_) \_\_\_\_\_  
**Confirmation Email:** \_\_\_\_\_

**Purchase Order #:** \_\_\_\_\_  
**Ship To:** \_\_\_\_\_  
**Attention:** \_\_\_\_\_  
**Address:** \_\_\_\_\_  
**Address:** \_\_\_\_\_  
**Ship To City:** \_\_\_\_\_ **ST:** \_\_\_\_\_ **ZIP:** \_\_\_\_\_  
**Ship To Phone:** \_\_\_\_\_

*The HCPCS CODES are not intended to be, nor should be considered billing or legal advice. They are provided based on code verification through PDAC or our interpretations of billing. The use of HCPCS CODES does not ensure coverage or payment. For coverage information, please verify the appropriate payer's coverage policy*

**S-ERGO ATX (KM-9000Series):**

*S-ERGO ATX Ultra lightweight Ergonomic Wheelchair std features: S-shape ergonomic seat frame, fixed footrest, quick release rear wheels, flip back armrest, composite side panels, tension adjustable seat and back sling, calf strap, frame weight 15.4 lb, max weight capacity 275 lb.*

BASE FRAME				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	S-ATX-1415	K0005	14x15 Foldable Ergonomic Active Wheelchair	\$2,877.00
<input type="checkbox"/>	S-ATX-1416	K0005	14x16 Foldable Ergonomic Active Wheelchair	\$2,877.00
<input type="checkbox"/>	S-ATX-1616	K0005	16x16 Foldable Ergonomic Active Wheelchair	\$2,877.00
<input type="checkbox"/>	S-ATX-1618	K0005	16x18 Foldable Ergonomic Active Wheelchair	\$2,877.00
<input type="checkbox"/>	S-ATX-1816	K0005	18x16 Foldable Ergonomic Active Wheelchair	\$2,877.00
<input type="checkbox"/>	S-ATX-1818	K0005	18x18 Foldable Ergonomic Active Wheelchair	\$2,877.00
FRAME COLOR				
<input type="checkbox"/>	LGCC		Lamborghini Verde Ithaca Lime Green	N/A
<input type="checkbox"/>	TOCC		Tango Orange	N/A
<input type="checkbox"/>	PWCC		Fashion White	std
<input type="checkbox"/>	SBCC		Diamond Black	std
SEAT CUSHIONS				
<input type="checkbox"/>	CU-ERGO-18	E2601	18x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-ERGO-16	E2601	16x17 seat memory foam 2" cushion	\$176.00
<input type="checkbox"/>	CU-FT-18	E2601	18x16 seat foam cushion	\$72.00
<input type="checkbox"/>	CU-FT-16	E2601	18x16 seat foam cushion	\$72.00

SEAT CUSHIONS cont.				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	CU-GFR-18	E2601	18x16 seat gel cushion	\$98.00
<input type="checkbox"/>	CU-GFR-16	E2601	18x16 seat gel cushion	\$98.00
<input type="checkbox"/>	US-SATX-1415		14x15 ½" Anti Bacterial seat cushion	std
<input type="checkbox"/>	US-SATX-1416		14x16 ½" Anti Bacterial seat cushion	std
<input type="checkbox"/>	US-SATX1616		16x16 ½" Anti Bacterial seat cushion	std
<input type="checkbox"/>	US-SATX1618		16x18 ½" Anti Bacterial seat cushion	std
<input type="checkbox"/>	US-SATX1816		18x16 ½" Anti Bacterial seat cushion	std
<input type="checkbox"/>	US-SATX1818		18x18 ½" Anti Bacterial seat cushion	std
BACK CUSHIONS				
<input type="checkbox"/>	CU-ERGO-18-back	E2611	18" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	CU-ERGO-16-back	E2611	16" wide memory foam 2" back cushion	\$116.00
<input type="checkbox"/>	BKF-2017	E2611	20" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-1817	E2611	18" wide back foam cushion	\$98.00
<input type="checkbox"/>	BKF-1617	E2611	16" wide back foam cushion	\$98.00
<input type="checkbox"/>	UB-SATX14		14" wide ½" Anti bacterial back cushion	std
<input type="checkbox"/>	UB-SATX16		16" wide ½" Anti bacterial back cushion	std
<input type="checkbox"/>	UB-SATX18		18" wide ½" Anti bacterial back cushion	std
REAR WHEELS				
<input type="checkbox"/>	RW-115		24x1 Pneumatic wheel with composite Ergonomic hand rim	Std
<input type="checkbox"/>	SPINLX		24x1 Spinergy LX 12 spoke, Pneumatic Kenda Tire and Tube w/ anodized polished handrim	\$895.00
<input type="checkbox"/>	SPINSPOX		24x1 Spinergy SPOX 18 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$695.00
<input type="checkbox"/>	SPINWIRE		24x1 Spinergy WIRE 30 spoke, Pneumatic Kenda Tire and Tube w/anodized polished handrim	\$295.00
KENDA TIRE COLOR/SPINERGY SPOKE COLOR (selection applies only to Spinergy Wheel)				
<input type="checkbox"/>	KENDAYLW		Yellow Kenda tire	Std
<input type="checkbox"/>	KENDARED		Red Kenda tire	Std
<input type="checkbox"/>	KENDABLU		Blue Kenda tire	Std
<input type="checkbox"/>	SPINSPYLW		Yellow Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPRED		Red Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPBLU		Blue Spinergy Spoke color	Std
<input type="checkbox"/>	SPINSPBLK		Black Spinergy Spoke color	Std
FRONT SEAT HEIGHT				
<input type="checkbox"/>	SHATX17		17" front seat height: 5x1 caster necessary	std
<input type="checkbox"/>	SHATX18		18" front seat height: 6x1 caster necessary	std
<input type="checkbox"/>	SHATX19		19" front seat height: 6x1 caster necessary	std
<input type="checkbox"/>	SHATX20		20" front seat height: 6x1 caster necessary	std

REAR SEAT HEIGHT				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	RHATX16		16" rear seat height	std
<input type="checkbox"/>	RHATX17		17" rear seat height	std
<input type="checkbox"/>	RHATX18		18" rear seat height	std
<input type="checkbox"/>	RHATX20		20" rear seat height	std
CENTER OF GRAVITY				
<input type="checkbox"/>	COGATXCOM		5 degrees rear wheel position: Stable users	std
<input type="checkbox"/>	COGATXSTD		0 degrees rear wheel position: Standard	std
<input type="checkbox"/>	COGATXACT		-5 degrees rear wheel position: Active users	std
BACK CANE HEIGHT				
<input type="checkbox"/>	BHATX105		10.5 inches back post height	std
<input type="checkbox"/>	BHATX11		11 inches back post height	std
<input type="checkbox"/>	BHATX12		12 inches back post height	std
<input type="checkbox"/>	BHATX13		13 inches back post height	std
<input type="checkbox"/>	BHATX14		14 inches back post height	std
<input type="checkbox"/>	BHATX15		15 inches back post height	std
BACK CANE ANGLE				
<input type="checkbox"/>	BAATX85		85 degrees back post angle: Active users	std
<input type="checkbox"/>	BAATX90		90 degrees back post angle: Standard	std
<input type="checkbox"/>	BAATX95		95 degrees back post angle: Comfortable	std
BACK CANE HANDLES				
<input type="checkbox"/>	BCHATX0		Omit push handles	std
<input type="checkbox"/>	BCHATX1		Push handle with non slip grip: Back cane height set to 14"	\$150.00
<input type="checkbox"/>	BCHATX2		Push handle with foldable non slip grip: Back cane height set to 14"	\$180.00
FOOTPLATE HEIGHT				
<input type="checkbox"/>	FHATX14		14" footplate height: seat sling to footplate	std
<input type="checkbox"/>	FHATX15		15" footplate height: seat sling to footplate	std
<input type="checkbox"/>	FHATX16		16" footplate height: seat sling to footplate	std
<input type="checkbox"/>	FHATX17		17" footplate height: seat sling to footplate	std
FRONT CASTERS (seat height dependent)				
<input type="checkbox"/>	C51LTALU		5x1 flat free polyurethane casters with polished aluminum hub	\$192.00
<input type="checkbox"/>	C51LT		5x1 flat free polyurethane casters composite hub	std
<input type="checkbox"/>	C61LT		6x1 flat free polyurethane casters composite hub	std
ARMREST ASSEMBLY				
<input type="checkbox"/>	AR9000FB		Flip back tubular armrest assembly with pad	std
<input type="checkbox"/>	AR9000NA		Omit armrest assembly	std

OPTIONS & ACCESSORIES				
✓	SKU#	HCPCS	Description	Price
<input type="checkbox"/>	SB22	E0978	Clamp style seat belt 2 piece	\$19.00
<input type="checkbox"/>	SB88	E0978	Wrap around seat belt 1 piece	\$19.00
<input type="checkbox"/>	SB99-48	E0978	Auto Style seat belt 2 piece	\$24.00
<input type="checkbox"/>	BT-115	E0961	Brake tip-Extending	Std
<input type="checkbox"/>	HR-FLD-115	E0955 + E1028	16-18" width Foldable Headrest with hardware	\$327.00
<input type="checkbox"/>	AT-9000	E0971	Anti Tippers-pair	\$80.00
<input type="checkbox"/>	LS16	K0038	16" width leg strap-calf support	\$11.00
<input type="checkbox"/>	LS18	K0038	18" width leg strap-calf support	\$11.00
<input type="checkbox"/>	LS20	K0038	20" width leg strap-calf support	\$11.00
<input type="checkbox"/>	PSG25F60-BLK	E1015	Frog Legs front fork suspension pair. Patient weight <160 lb	\$299.00
<input type="checkbox"/>	PSG25F61-BLK	E1015	Frog Legs front fork suspension pair. Patient weight >160 lb	\$299.00
<input type="checkbox"/>	CP2	E1399	Universal carry pouch-small	\$11.00
<input type="checkbox"/>	CP4	E1399	Universal carry pouch-large	\$19.00
<input type="checkbox"/>	CANO-115		UV protectant Canopy for 7/8" handle frame	\$248.00
<input type="checkbox"/>	PB-115-16		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	PB-115-18		Foldable rigidifying push bar for 16" width	\$88.00
<input type="checkbox"/>	PB-115-20		Foldable rigidifying push bar for 16" width	\$88.00