



S-Ergo 305

Top End Lightweight Ergonomic Manual Wheelchair

- Medicare code K0005
- High strength, weighs only 29 lbs. (w/o footrests)
- Only 22lbs (w/ wheels & footrests removed)
- Anti-Bacterial Upholstery/ Cushion
- Folding backrest, for transporting chair in car, bus, trips, ect.
- Flip back armrests w/ wider contoured arm pads for maximum comfort
- Swing in & away footrest for maximum safety while entering or exiting the chair
- “Tube-in Center” foot-plate, assures better side leg support
- 24” Quick Release rear spoke/polyurethane/high profile/ flat free wheels & 7x1” front casters
- 16x17 or 18x17 Seat width
- Optional: Height adjustable seat per request
- Frame Color: Pearl Silver & Rose Red
- Upholstery: Silver/Black mesh bottom/ AEGIS top
- Weight Capacity of 250 lbs

Features for S-Ergonomic Series

- Ergonomic, ultra lightweight specialized frame
- Pressure relief, stabilization, and anti-slipage seating system
- Upholstery is equipped with AEGIS anti-biotic technology, odor/stain/deterioration control.
- Wider/Concaved arm pad, providing comfort & support

Model No.	seat width (in.)	seat depth (in.)	armrest height (in.)	seat height (in.)	back height (in.)	overall height (in.)	overall open width (in.)	folded width (in.)	overall length (in.)	weight w/o riggings (lb.)	weight capacity (lb.)
S-305-16	16	17	9-11	19	17	38	24	12	43	29	250
S-305-18	18	17	9-11	19	17	38	26	12	43	29	250

**Available Options ~ Anti-Tipper, Seat belt, Extra Upholstery, Memory Foam Cushion, Ergonomic/Contoured Cushion Reflecting Spoke Guard, Frog Leg Suspension, Magnesium Wheel, Vulcan Grip Wheel, Angle/Height/Length Adjustable Headrest*

Folded Dimensions: 33” x 12” x 28” (L x W x H)



KARMAN HEALTHCARE Inc.

19255 San Jose Avenue. City of Industry, CA 91748-1418
 Order. 1-800-805-2762 Tel. 626-581-2235 Fax. 626-581-2335